



Cohesive Bandages

INSTRUCTIONS FOR USE

INTENDED USE / PURPOSE

Designed to provide controlled compression and support for soft tissue injuries, secure dressings, and stabilize body parts following trauma or medical procedures:

- Compression for sprains, strains, and mild swelling
- Securing dressings, splints, or tubing
- Support for joints or muscles during recovery

CONTRAINDICATIONS

Do not use directly on open wounds without a primary dressing. Avoid use if the patient has poor circulation or if skin irritation occurs.

SUGGESTED APPLICATION BY SIZE

- 1 in (2.5 cm): Fingers, toes; securing dressings, catheters, or tubing
- 2 in (5 cm): Fingers, wrists, or small dressings
- 3 in (7.5 cm): Hands, forearms, ankles
- 4 in (10 cm): Elbows, knees, or smaller limbs
- 6 in (15 cm): Thighs, large limbs, or torso areas

DIRECTIONS FOR USE

1. Prepare the area: Ensure skin is clean and dry. Apply a sterile dressing if needed.
2. Start wrapping: Anchor the bandage with one to two turns, overlapping each layer by 50%.
3. Adjust tension: Apply gentle, even pressure—tight enough for support but not to restrict circulation.
4. Secure: Press the final layer lightly so it adheres to itself. No clips or tape required.
5. Check circulation: Ensure the area beyond the bandage is warm and has normal color and sensation.

WARNINGS AND PRECAUTIONS

- Do not stretch excessively during application.
- Reapply if the bandage loosens or becomes wet.
- For single patient use only.

STORAGE

Store in a cool, dry place away from direct sunlight and moisture.

CONCORDANCE
Sourcing Solutions

concordancehealthcare.com

Manufactured for and distributed by:
Concordance Healthcare Solutions, LLC. Tiffin, OH 44883
1-855-241-1773

CALIFORNIA PROPOSITION

⚠ 65 WARNING ⚠

WARNING: This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

For more information:
www.P65Warnings.ca.gov