



Elastic Bandages

INSTRUCTIONS FOR USE



Not made with natural rubber latex
No está hecho de látex de caucho natural



Single Use
Uso único

INTENDED USE / PURPOSE

Used to support and compress soft tissue injuries, reduce swelling, and secure dressings. Suitable for use in both first aid and clinical care to manage strains, sprains, or edema.

- Support for muscle or joint injuries
- Compression therapy for swelling or varicose veins
- Post-surgical or trauma support
- Securing dressings or splints

CONTRAINDICATIONS

Do not apply directly over open wounds without a protective dressing. Avoid if the patient has impaired circulation, severe arterial disease, or known sensitivities to elastic materials.

SUGGESTED APPLICATION BY SIZE

- 2 in (5 cm): Fingers, hands, wrists
- 3 in (7.5 cm): Forearms, ankles
- 4 in (10 cm): Knees, elbows, or smaller limbs
- 6 in (15 cm): Thighs, large limbs, or torso areas

DIRECTIONS FOR USE

1. Prepare the area: Clean and dry the skin; apply a sterile dressing if required.
2. Anchor: Begin with one to two turns to anchor the bandage.
3. Wrap: Overlap each layer by about 50%. Maintain consistent, moderate tension.
4. Secure: Use provided clips or fasteners to hold the end in place.
5. Check fit: The bandage should be snug but not tight—ensure circulation remains normal.

WARNINGS AND PRECAUTIONS

- Do not apply over tight or bulky dressings.
- Remove immediately if numbness, pain, or discoloration occurs.
- Wash and air-dry between uses, following manufacturer's instructions.
- Single-patient use recommended for infection control.

STORAGE

Store in a dry environment at room temperature. Avoid direct sunlight and heat sources.

CONCORDANCE
Sourcing Solutions

concordancehealthcare.com

Manufactured for and distributed by:
Concordance Healthcare Solutions, LLC. Tiffin, OH 44883
1-855-241-1773

CALIFORNIA PROPOSITION

⚠ 65 WARNING ⚠

WARNING: This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

For more information:
www.P65Warnings.ca.gov