

INTRODUCING



CRYOSURGERY WITHOUT THE GUESSWORK

- Nitrous Oxide OR Carbon Dioxide Gas—First and only cryo product that offers the flexibility to use either cryogen based on your patients' needs.
- * Advanced Treatment Wand—Now you can treat your patients in any direction, at any angle, on any part of the body.
- * Preset Dispense Times—Choose between 5-second or 10-second treatment times and minimize the risk of under- or over-treating.
- ** Gas Level Indicator—Our new lighted indicator ensures you will know when to order a new tank.






CRYOSURGERY JUST GOT A WHOLE LOT SMARTER!

For the first time ever, you can choose between two cryogens—Nitrous Oxide (-89°C) or Carbon Dioxide (-79°C)—to get great first-time outcomes and minimize damage to surrounding skin.

N₂O is an ideal cryogen for most common lesions such as verruca vulgaris, angiomas, seborrheic keratosis and more.*

What cryogen is best for cosmetic lesions such as skin tags, age and sun spots on the hands and face? The answer is CO₂ which provides a more superficial treatment and minimizes your patients' downtime and recovery.

FREEZING POINT COMPARISON



Item No.	Product Code	CryoLab® Medical Product Description
200-1000	CryoLab Medical Device Unit	Kit Contains: 1 Canister of cryogen 1 N ₂ O OR CO ₂ 1 12V/20W Power Supply 1 90" Power Cord 1 Quick Start Guide 1 Instructions for Use One-year Limited Warranty*
000-2730	Replacement N ₂ O Cylinder	Contains: 1 Canister of N ₂ O cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label
000-2731	Replacement CO ₂ Cylinder	Contains: 1 Canister of CO ₂ cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label



Contact your Prime Vendor TODAY!