



## Product Discontinuation Notice – ACTION REQUIRED

RE: Discontinuation of Select Laryngoscope Replacement Lamps

December 10, 2025

Dear Valued Customer,

As part of our business process, AirLife routinely evaluates its product portfolio and occasionally discontinues the manufacturing and distribution of certain products and materials. AirLife acknowledges the obligation to promptly inform our customers of the intent to exit products to enable our customers to prepare for alternatives and assist in their product portfolio renewal process. This letter serves as a formal notification that AirLife has discontinued production of the product(s) found on ATTACHMENT A.

**Carefully review this information and direct this notice to your procurement team or item discontinuation team to remove these items from your systems and transition to the suggested replacements.**

We value your relationship, and our intent is to provide ample notice to mitigate the inconvenience and impact this may cause. Our Sales Representatives are available to help you with the recommended substitute options and to answer any questions. You can also reach out to our Customer Service Team at [info@myAirLife.com](mailto:info@myAirLife.com) or 800.433.2797.

Sincerely,

*Nick Soler*

Nick Soler  
Associate Director of Product Management

ENC:  
Attachment A, Attachment B



ATTACHMENT A

CURRENT VENDOR MATERIAL				REPLACEMENT MATERIAL		
Product SKU Number	Product Description	Sales Discontinuation Notice Date	Current Inventory	Substitute SKU Number	Substitute Description	Additional Comments
5-0238-60	F/O Conventional Lg Thread	Dec 5, 2025	31 EA	N/A	N/A	No Substitute Available
5-0234-14	Conv Reflector Lamp Large Pk/5	Dec 5, 2025	0 EA	N/A	N/A	No Substitute Available
5-0240-59	F/O Lamp G-L Black Base 3.5mm	Dec 5, 2025	0 EA	5-0240-61 5-0240-62 5-0240-63	See Below	Substitute lamps compatible with specific handles
5-0240-52	F/O Lamp Red Base 10/32	Dec 5, 2025	11 EA	5-0240-61 5-0240-62 5-0240-63	See Below	Substitute lamps compatible with specific handles

The manufacturing date of your fiber-optic laryngoscope handles can be identified by the two-letter stamp on the head of the laryngoscope. The first letter is upper case and identifies the year while the second letter is lower case and identifies the month. We are no longer able to supply replacement lamps for handles with upper case letters prior to M.

SKU 5-0240-59 is compatible with fiber-optic handles with a lot number that started with capital prior to H.

SKU 5-0240-52 is compatible with fiber-optic handles with a lot number that started with capital H through L.



## ATTACHMENT B

The chart below outlines which replacement lamps are compatible with our current line of fiber-optic laryngoscope handles made in 2012 and later. This is for all lots stamped with a first letter of a capital M or later:

	<b>SKU:</b>		5-0240-61	5-0240-62	5-0240-63
	<b>Cartridge Color:</b>		Blue	Orange	Yellow
	<b>Material</b>	<b>Size</b>	Halogen Xenon	Econo LED	2W LED
5-0236-09	Chrome	Medium	X		
5-0236-10		Penlite	X		
5-0236-11		Stubby	X		
5-0236-15		Large	X		
5-0236-30	Stainless Steel	Penlite	X		
5-0236-31		Stubby	X		
5-0236-35		Large	X		
5-0236-36		Micro	X		
5-0236-38		Mini	X		
5-0236-39		Medium	X		
5-0236-90	Chrome	Penlite		X	
5-0236-91		Stubby		X	
5-0236-99		Medium		X	
5-0236-80	Chrome	Penlite			X
5-0236-89		Medium			X
5-0246-80	Stainless Steel	Penlite			X
5-0246-81		Stubby			X
5-0246-86		Micro			X
5-0246-88		Mini			X
5-0246-89		Medium			X