

September 2025

PRODUCT MODIFICATION

VACUETTE® BLOOD CULTURE HOLDER

Dear Valued Customer,

We would like to inform you of an exciting upcoming product change. At Greiner Bio-One, it is essential for us to keep our products up to date, compliant, and aligned with current market needs. We would like to inform you about a recent design review of the VACUETTE® Blood Culture Holder. The following changes have been made to the design:

Addition of prominent single-use symbol:

This visual indication ensures that the holder is not intended for re-use, which is crucial to prevent contamination from potentially infected blood samples.

VACUETTE® logo replaced by Greiner Bio-One logo:

This change allows for immediate visual recognition of Greiner Bio-One manufactured products, which when used as part of a system provide optimal compatibility and functioning.



New indication for
single use only



Previous Blood Culture Holder (left) and
new Blood Culture Holder with Greiner
Bio-One logo (right)

Furthermore, a minor change has been made to the holder thread to align all Greiner Bio-One holder products, optimizing product compatibility across the VACUETTE® product line.

The following items are affected:

Item No.	Description	Qty. (bag/case)
450181	VACUETTE® Blood Culture Holder	40 / 320

Greiner Bio-One North America

4238 Capital Drive, Monroe, North Carolina 28110


PREANALYTICS PHONE 1-888-286-3883 | **BIOSCIENCE PHONE** 1-800-884-4703

FAX 1-800-726-0052 | **E-MAIL** office@us.gbo.com | **WEBSITE** www.gbo.com

This change will apply to all products manufactured from October 2025 onwards. However, it is possible that you may still receive inventory from earlier production lots manufactured prior to the change.

As always, the Greiner Bio-One Technical Services team is available to support you in the use of all our products at (800) 515-8112 or patech@gbo.com. You may also contact Customer Service at (888) 286-3883 or your local Greiner Bio-One Account Manager for assistance.

Best regards from Greiner Bio-One,



Meredith Millay
Product Manager - Preanalytics

Greiner Bio-One North America

4238 Capital Drive, Monroe, North Carolina 28110

PREANALYTICS PHONE 1-888-286-3883 | **BIOSCIENCE PHONE** 1-800-884-4703

FAX 1-800-726-0052 | **E-MAIL** office@us.gbo.com | **WEBSITE** www.gbo.com