

New Product Announcement Biogel® PI Micro Indicator® System



Dear Valued Distributor:

April 12, 2024

Mölnlycke Health Care, manufacturer of Biogel surgical gloves, is please to introduce a new synthetic double-gloving system, the **Biogel PI Micro Indicator System** that offers safety, comfort, efficiency and sustainability – all in one convenient package.

The Biogel® PI Micro Indicator® System offers:

Improved convenience

- A single pack containing two pairs of gloves designed and sized specially to glide on easily and work perfectly together.

Improved sustainability


- Convenient double-packed gloves with 50% less plastic packaging by eliminating a sterile pouch and packing in smaller boxes

Improved protection

- Consists of a blue, polyisoprene surgical indicator underglove and a straw-colored overglove, creating a patented puncture indication system proven to provide clear, fast and large¹ visual alert to the wearer reducing the risk of occupational exposure

Improved comfort

- Biogel® PI Micro Indicator® System gloves are extra thin and conform to the hand for maximum tactile sensitivity when double gloving.²
- ½ size larger underglove enhances comfort



Biogel® PI Micro Indicator® System (48355–48385) features:

- Sterile, synthetic, polyisoprene surgical glove (manufactured without natural rubber latex)
- Contains blue underglove ½ size larger than overglove, tested for use with chemotherapy agents
- Contains straw color overglove, tested for use with chemotherapy agents
- Made with a curved former for a more natural fit
- Formulated for improved sensitivity even when double gloving
- Standard grip and anti-slip overglove cuff

Recommended use

- Recommended for all surgical procedures and especially those when extra tactile sensitivity is desired even when double gloving.

Please contact your local representative or call the Mölnlycke Customer Service Center at 1-800-843-8497 for more information. Thank you for your support and we look forward to continuing to help meet your surgical needs.

Regards,
Corinne Schmid

A handwritten signature in black ink that reads 'Corinne Schmid'.

Senior Marketing Director, Biogel®
Mölnlycke

¹Summary of Indication Performance of Biogel Indicator Systems versus Competitors' Double Gloving Combinations. Mölnlycke Health Care, 2020. Data on file

²Summary of Technical Documents. Mölnlycke Health Care. Data on File.



Ordering Information - Biogel® PI Micro Indicator® System (Ordering UOM is by CASE)

Re-order #	Size	Pairs per pouch	Pairs / Box	# Pairs / Case	List Price / Case \$
48355	5.5	2	2 x 25 (50)	2x100 (200)	\$504
48360	6	2	2 x 25 (50)	2x100 (200)	\$504
48365	6.5	2	2 x 25 (50)	2x100 (200)	\$504
48370	7	2	2 x 25 (50)	2x100 (200)	\$504
48375	7.5	2	2 x 25 (50)	2x100 (200)	\$504
48380	8	2	2 x 25 (50)	2x100 (200)	\$504
48385	8.5	2	2 x 25 (50)	2x100 (200)	\$504

ORDER #	UOM	# of PCS	EAN/UPC BARCODE	GTIN	Material Content	Sterility
48355	PR	2	7323190139255	7323190139255	Polyisoprene	Sterile
48355	BX	2 x 25	7323190091188	7323190091188	Polyisoprene	Sterile
48355	CS	2 x 100	7323190091171	7323190091171	Polyisoprene	Sterile

48360	PR	2	7323190139262	7323190139262	Polyisoprene	Sterile
48360	BX	2 x 25	7323190091225	7323190091225	Polyisoprene	Sterile
48360	CS	2 x 100	7323190091218	7323190091218	Polyisoprene	Sterile

48365	PR	2	7323190139279	7323190139279	Polyisoprene	Sterile
48365	BX	2 x 25	7323190091263	7323190091263	Polyisoprene	Sterile
48365	CS	2 x 100	7323190091256	7323190091256	Polyisoprene	Sterile

48370	PR	2	7323190139286	7323190139286	Polyisoprene	Sterile
48370	BX	2 x 25	7323190091300	7323190091300	Polyisoprene	Sterile
48370	CS	2 x 100	7323190091294	7323190091294	Polyisoprene	Sterile

48375	PR	2	7323190139293	7323190139293	Polyisoprene	Sterile
48375	BX	2 x 25	7323190091348	7323190091348	Polyisoprene	Sterile
48375	CS	2 x 100	7323190091331	7323190091331	Polyisoprene	Sterile

48380	PR	2	7323190139309	7323190139309	Polyisoprene	Sterile
48380	BX	2 x 25	7323190091386	7323190091386	Polyisoprene	Sterile
48380	CS	2 x 100	7323190091379	7323190091379	Polyisoprene	Sterile

48385	PR	2	7323190139316	7323190139316	Polyisoprene	Sterile
48385	BX	2 x 25	7323190091461	7323190091461	Polyisoprene	Sterile
48385	CS	2 x 100	7323190091454	7323190091454	Polyisoprene	Sterile



Storage Information – Biogel® PI Micro Indicator® System (Ordering UOM is by CASE)

ORDER #	UOM	LENGTH	WIDTH	HEIGHT	DIMENSION UNIT	VOLUME	VOLUME UNIT	WEIGHT	WEIGHT UNIT
48355	PR	11.34	5.91	-	IN	-	CU IN	1.78	OZ
48355	BX	11.34	5.98	7.91	IN	0.310	CU FT	3.19	LB
48355	CS	23.15	12.60	8.66	IN	1.459	CU FT	14.48	LB

48360	PR	11.34	5.91	-	IN	-	CU IN	1.88	OZ
48360	BX	11.34	5.98	7.91	IN	0.310	CU FT	3.33	LB
48360	CS	23.15	12.60	8.66	IN	1.459	CU FT	15.07	LB

48365	PR	11.34	5.91	-	IN	-	CU IN	1.96	OZ
48365	BX	11.34	5.98	7.91	IN	0.310	CU FT	3.46	LB
48365	CS	23.15	12.60	8.66	IN	1.459	CU FT	15.58	LB

48370	PR	11.34	5.91	-	IN	-	CU IN	2.11	OZ
48370	BX	11.34	5.98	7.91	IN	0.310	CU FT	3.70	LB
48370	CS	23.15	12.60	8.66	IN	1.459	CU FT	16.54	LB

48375	PR	11.34	5.91	-	IN	-	CU IN	2.24	OZ
48375	BX	11.34	5.98	7.91	IN	0.310	CU FT	3.90	LB
48375	CS	23.15	12.60	8.66	IN	1.459	CU FT	17.35	LB

48380	PR	11.34	5.91	-	IN	-	CU IN	2.35	OZ
48380	BX	11.34	5.98	7.91	IN	0.310	CU FT	4.07	LB
48380	CS	23.15	12.60	8.66	IN	1.459	CU FT	18.01	LB

48385	PR	11.34	5.91	-	IN	-	CU IN	2.48	OZ
48385	BX	11.34	5.98	7.91	IN	0.310	CU FT	4.27	LB
48385	CS	23.15	12.60	8.66	IN	1.459	CU FT	18.82	LB

Storage Instructions: Store in a dry place at a temperature of 5-25°C, away from sources of heat or direct sunlight.