



INTRODUCING OUR NEW PRODUCT: **Linearity FD Urine Chemistry** **for Beckman AU**



Linearity FD Urine Chemistry for Beckman AU (Order Number: K829M-10) is intended to simulate human patient samples for determining linearity, calibration verification, and the verification of reportable range for its respective analytes.

Package Size: 10 x 3 mL

Type: Freeze Dried

Levels: 10

Open Vial Stability: 5 days when stored at 2-8°C


Analytes: Set 1: Phosphorous, Urea Nitrogen, Uric Acid

Set 2: Amylase, Calcium, Chloride, Creatinine, Glucose (Urine), Glucose (CSF), Magnesium, Microalbumin (Urine), Microalbumin (CSF), Potassium, Protein, Sodium

FREE DATA REDUCTION SERVICES

AUDITOR QC

- Free data reduction software
- Peer group data
- Inspector ready charts and graphs
- Real-time online data analysis



AUDIT FAXBACK

- Fax in your calibration verification / linearity results and Audit will fax back reports.