

DISPOSABLE PULSE OXIMETRY SENSORS


The internal surfaces of pulse oximetry sensors are frequently soiled with various materials, including microorganisms. Since reusable pulse oximetry sensors come into direct contact with a patient's skin, and are used with multiple patients, they are a potential source of cross contamination and HAI transmission.

At Curbell, we understand that hospital-acquired infections (HAIs) are among the top ten leading causes of death in the United States, which is why we're pleased to offer a line of disposable pulse oximetry sensors. Utilizing these sensors is a simple way to help reduce complex cross-contamination issues and enhance patient care.

- ▶ Comfortable, form-fitting sensors
- ▶ Designed for single patient use to provide a solution to your facility's infection control program
- ▶ Available in adult, neonate, pediatric, and infant
- ▶ Nellcor Oximax™, Non-Oximax, and Nihon Kohden compatibility
(Oximax connectors are shown here — see the chart below for all available styles)



Adhesive Film



Adhesive Fabric



Non-Adhesive Flocking

SIZE/WEIGHT	TECHNOLOGY	CROSS REFERENCE PN	LENGTH	MATERIAL	CASE QTY	CURBELL PN
Adult (>30 kg)	Nellcor Non-OxiMax	D-25	39"	Adhesive Film	24	28514
	Nellcor OxiMax	MAXA and MAXA-NS	20"	Adhesive Film	24	28515
		MAXAL	39"	Adhesive Film	30	28157
	Nihon Kohden	TL-271T	39"	Adhesive Film	24	28534
		TL-271T3	63"	Adhesive Film	24	28172
Infant (3-20 kg)	Nellcor Non-OxiMax	I-20	39"	Adhesive Fabric	24	28517
	Nellcor OxiMax	MAXI and MAXI-NS	20"	Adhesive Fabric	24	28518
		MAXIL	39"	Adhesive Fabric	24	28519
	Nihon Kohden	TL-274T3	63"	Adhesive Fabric	24	28174
Neonatal/Adult (<3 kg or >40 kg)	Nellcor Non-OxiMax	N-25	39"	Adhesive Fabric	24	28523
		SC-NEO	39"	Non-Adhesive Flocking	30	26239
	Nellcor OxiMax	MAXN and MAXN-NS	20"	Adhesive Fabric	24	28524
		MAXNL	39"	Adhesive Fabric	24	28525
		SC-A	39"	Non-Adhesive Flocking	30	27092
		SC-NEOL	39"	Non-Adhesive Flocking	24	28522
	Nihon Kohden	TL-273T	39"	Adhesive Fabric	24	28535
		TL-273T3	63"	Adhesive Fabric	24	28180
Pediatric (10-50 kg)	Nellcor Non-OxiMax	D-20	39"	Adhesive Film	24	28531
	Nellcor OxiMax	MAXP	39"	Adhesive Film	24	28532
		MAXPL	39"	Adhesive Film	24	28533
	Nihon Kohden	TL-272T	39"	Adhesive Film	24	28536
		TL-272T3	63"	Adhesive Film	24	28178

NOTE: All Curbell disposable SPO2 sensors are non-sterile. Nellcor and OxiMax are trademarks of Nellcor, Inc. Nihon Kohden is a trademark of Nihon Kohden Corporation.

©2025 Curbell Medical Products, Inc. MAP1795E