

SKU	Description	Safety Version	Differences	Cardinal Sub	BD Sub
SYR100167	16G x 1.5"	No safety sub	same needle + safety mechanism	None	B-D305198
SYR100185	18G x 1"	SSN100185	same needle + safety mechanism	SWD250040	B-D305195
SYR100187	18G x 1.5"	SSN100187	same needle + safety mechanism	SWD818112	B-D305196
SYR100195	19G x 1"	SSN100195	same needle + safety mechanism	SWD819100	None
SYR100197	19G x 1.5"	No safety sub	same needle + safety mechanism	SWD819112	None
SYR100205	20G x 1"	SSN100205	same needle + safety mechanism	SWD250123	None
SYR100207	20G x 1.5"	SSN100207	same needle + safety mechanism	SWD820112	None
SYR100215	21G x 1"	SSN100215	same needle + safety mechanism	SWD250172	B-D305165
SYR100225	22G x 1"	SSN100225	same needle + safety mechanism	SWD822100	None
SYR100227	22G x 1.5"	SSN100227	same needle + safety mechanism	SWD822112	None
SYR100235	23G x 1"	SSN100235	same needle + safety mechanism	None	B-D305145
SYR100253	25G x 5/8"	SSN100253	same needle + safety mechanism	SWD250313	B-D305122
SYR100255	25G x 1"	SSN100255	same needle + safety mechanism	None	B-D305125
SYR100257	25G x 1.5"	SSN100257	same needle + safety mechanism	SWD825112	B-D305127
SYR100302	30G x 1"	SSN100302	same needle + safety mechanism	SWD830340	B-D305106

