

Acute Care & Monitoring

6135 Gunbarrel Ave

Boulder, CO 80301

www.medtronic.com



RE: Mon-a-therm™ Foley catheter with temperature sensor 400TM item code change

Dear Valued Customer,

Thank you for your trust and continued use of Medtronic products. At Medtronic, we value continuous improvement and strive to deliver sustainable technologies that enhance patient care.

We are pleased to announce the release and immediate availability of a new version of the Mon-a-therm™ Foley catheter with temperature sensor 400TM, which will replace the older version. The Mon-a-therm™ Foley catheter is designed to provide urinary bladder drainage while simultaneously measuring a patient's core body temperature at the insulated site.

The following are comparative details between the two versions of the product. It is important to note that the function of the previous catheter and the new replacement remain the same. Clinicians might notice only minor cosmetic differences, such as tubing color and packaging.

Old Mon-a-therm™ Foley catheter codes		New Mon-a-therm™ Foley catheter codes	
			
Code	Description	Code	Description
90051	16FR Foley Catheter with Temp	90051T	16FR Foley Catheter with Temp
90052	18FR Foley Catheter with Temp	90052T	18FR Foley Catheter with Temp

90053	8FR Foley Catheter with Temp	90053T	8FR Foley Catheter with Temp
90054	10FR Foley Catheter with Temp	90054T	10FR Foley Catheter with Temp
90055	12FR Foley Catheter with Temp	90055T	12FR Foley Catheter with Temp
90056	14FR Foley Catheter with Temp	90056T	14FR Foley Catheter with Temp

Our dedicated team is here to support you. For any queries, clarifications, or more information, please do not hesitate to contact our customer support team at 1-800-722-8772 or your local Medtronic Acute Care and Monitoring field sales representative.

Sincerely,

Medtronic US Marketing
Acute Care & Monitoring